

# ExPendBox



SMALL IN SIZE.  
BIG ON CAPABILITY.



## ExpendoBox

Functionality and compact design make it easily integrable in any clinic.


Gloves, masks, paper, towels ... ExpendoBox is a dispenser for the disposable material that you usually need in your podiatry clinic.

Extremely functional and compact, it allows you to always have all the essential consumables for your daily work.



# PRACTICAL AND FUNCTIONAL WITH AN ATTRACTIVE DESIGN

# EXPENDBOX



The door has a folding system with brake for greater comfort.



Made of compact laminate with a white lacquered steel frame.



You can choose your dispenser with the details painted in the color of your chair.




Save space on your unit desktop thanks to the Expendbox storage capacity.





### Sharps container

Safely dispose of scalpels, needles, or other sharp or sharp materials.



### Gloves

It has two dispensers, one on the side and the other on the bottom to be able to dispense two types or two sizes of gloves.



### Face Masks

Mask dispenser for easy and hygienic dispensing.




### Paper hand roll

Continuous wick roll paper dispenser. Practical and hygienic.




### Towels

Viewing window for hygienic individual dispensing of sizes.



# EXPENDBOX

**Expendbox can contain:**

- 50 towels
- 2 glove boxes
- 1 box of masks
- 1 roll of hand drying paper
- 1 disposable scalpel container

**Closed measures:**

480mm x 360mm x 297mm

**Interior measures:**

467mm x 347mm x 267mm

**Weight:**

7 Kgs.